## Increase titer and improve product quality using perfusion culture



Pingfang Chen, Shan Liu, Rubin Jiang GenScript ProBio Inc., NJ USA 008854

## Abstract

Issues like low yield, aggregation, fragmentation, charge variant and hard to purify impurities are frequently found in upstream process for different types of biologics. Perfusion process can promote cell growth and improve product quality. In perfusion process, fresh media is pumped into bioreactor and product-containing media is exchanged out. In GenScript Probio platform, two modes of perfusion culture were established to better produce unstable molecules and molecules associates with hard to purify impurities.

Items	Static perfusion	Dynamic perfusion
Cell bleeding	YES	NO
Loss of cell and product yield	YES	NO
Operation	More complicated	Less complicated
Culture duration	Long (30 – 100 days)	Short (18 – 30 days)
Cell viability	Maintained well through out the culture	Decreased in the late stage of culture

## Static perfusion ■D28 ■D32 ■D36 ■D57 ■D61 ■D65 ■D69 VCD (10<sup>6</sup> cells/n VIA ■ D73 Percentage VVD 40 2 VVD 28 32 36 40 44 48 52 56 60 64 68 72 76 80 84 30 40 50 60 Time (day) Time (day) Perfusion showed higher main peak of Culture duration was up to 84 days. Productivity was about 28 times of The VCD and viability were maintained traditional FB (4 g/L vs 110 g/L). iCIEF and lower acidic species. well at the static phase. The quality attributes were stable during the static phase.

