# The Platform Strategy of Antibody Disulfide **Bond Reduction to Reduce the Risk of Antibody Process Development**




Shuo Tang, Jiaoli Tao, Qianqian Zhu, Ying Li GenScript ProBio Inc., NJ USA 008854

## ABSTRACT

The reduction of disulfide bonds poses a highly challenging issue in antibody production, as it can lead to a decrease in product purity and can impact the safety and efficacy of drugs. Therefore, it is necessary to develop and optimize process strategies to minimize disulfide bond reduction and ensure the quality, efficacy, and safety of the product. In recent years, with the emergence of more complex molecules such as bispecific and trispecific antibodies, the problem of molecular heterogeneity caused by incomplete or mispaired interchain disulfide bonds has become even more daunting.

Recently, ProBio has developed the ProBox<sup>™</sup> platform, which offers efficient solutions for antibody reduction. This strategy involves initially assessing the risk of molecular reduction and using the resulting score to guide the implementation of additional mitigation strategies. In summary, we present a risk assessment strategy for antibody disulfide bonds reduction, and the strategy to mitigate the impact of antibody reduction on subsequent purification processes.

## ANTIBODY DISULFIDE BOND REDUCTION PRINCIPLE, RISK ASSESSMENT AND SOLUTION STRATEGY




#### Ren, Tingwei, et al. Biotechnology and Bioengineering 118.8 (2021): 2829-2844.


First Reaction in Glycolysis

Antibody disulfide bond reduction is essentially an oxidation-reduction (redox) reaction that involves redox enzyme. Glutathione (GSH) and thioredoxin (Trx) system (comprising Trx, thioredoxin reductase [TrxR] and nicotinamide adenine dinucleotide phosphate (NADPH) are the known enzymes, and enzyme systems, that contribute to disulfide reduction. NADPH is generated from the pentose phosphate pathway and serves as an electron source in the disulfide bond reduction. Electrons first transfer from NADPH to TrxR and reduce the TrxR disulfide bond, then move to the oxidized Trx to form the reduced Trx, and finally reduce the disulfide bond. GSH catalyzes the disulfide bond reduction in a similar way as Trx system.

	Impact factor	Grading standards	score	percentage	weighted score (score*percentage)
		high sensitivity	10		
Reduction	Reducing agent test <sup>*1</sup>	moderately	-	20.9/	



## **CASE STUDY: FUSION PROTEIN**



Solutions

CE CDC ND /0/

Solutions:		Sample	CE-SDS-NR (%)			30			
			pre-peaks	Main Peak	post-peaks	Contraction of the second s			
Increase dissolved oxygen to consume	0	Fresh supernatant	2.5	88.8	8.8	and the second s			
NADPH	1	Place at RT for 7 h	36.0	56.7	7.3		0		
A distance to be a set of the ball that There is a the life of	2	Continuous air sparging 20 h	1.4	90.6	8.1				
Add metal ions to inhibit Trx activity	3	Insulation 24 h	40.9	56.4	2.7				
✓ Add oxidant consuming Trx	4	Add 1 mM CuSO <sub>4</sub> 27 h	2.1	89.7	8.2	Continuous			
	5	Add 1 mM cystamine 30 h	1.5	91.5	7.0	Continuous Air sparging 20 h	Insulate the air 24 h	Add 1 mM	Add 1 mM
							LITE all 24 II	CuSO₄ 27 h	cystamine 30 h

## CONCLUSION

We have established a comprehensive platform for antibody disulfide bond reduction, which assesses the risk of disulfide bond reduction during the early stages of antibody development. Based on the severity of the risk, we recommend appropriate solution strategies to mitigate potential risks in later stages of development. This platform has been successfully applied to antibody projects, ensuring smooth process development and obtaining highly stable therapeutic products.



www.genscriptprobio.com